

FoxP3 intracellular staining on microtiter plates

eBioscience Kit

1. We put the cell suspension on microtiter plate by 1×10^6 cells/well.
2. 1x wash with 200 μ l PBS (400g, 2min., 4°C)
3. Stain surface antigens (see Surface staining)
4. 1x wash with 200 μ l PBS
5. shake
6. resuspend cells in 50 μ l fresh Fixation/Permeabilization buffer:

Preparing Fix./Perm, buffer:

1 unit of Fixation/Permeabilization Concentrate + 3 units of Fixation/Permeabilization Diluent. (e.g. 100 μ l +300 μ l)

7. incubate 30min. 4°C
1x wash with 100 μ l 1x permeabilization buffer

Preparing 1x permeabilization buffer:

*1 unit of 10x Permeabilization Buffer + 9 units of dd H₂O
(e.g. 1ml 10x Perm. Buf. +9ml ddH₂O)*

8. Block cells with 50 μ l 2% rat serum diluted in 1x perm. buffer
9. Incubate 20min 4°C
10. Centrifuge, aspirate liquid
11. Dilute FoxP3 antibody at given titer in 1x perm. buffer
12. Stain with diluted FoxP3 antibody by 10 μ l per well
13. Incubate 30min., 4°C, in dark
14. 2x wash with 200 μ l 1x perm. buffer
15. Resuspend cells in 100 μ l 1x perm. buffer
16. analyse